

# EDSflow™

## Continuous Flow Effluent Decontamination System



### The Ultra-efficient Continuous Flow BioWaste Effluent Decontamination System

#### Overview

The Suncombe EDSflow™ Continuous Flow BioWaste Decontamination Kill System has been developed to provide an Ultra-efficient continuous decontamination system.

Following numerous clients' requests to develop a continuous based biowaste decontamination system that could be validated, using Suncombe tremendous experience of Biowaste decontamination systems and critical processing systems, we have leveraged our knowledge of batch and hybrid biowaste decontamination systems.

The EDSflow™ is a pre-designed system which treats effluent continuously. It consists of a buffer tank which buffers the effluent, ensuring that there is a controlled supply to the continuous system regardless of flow peaks and troughs.

Using steam, super-heated water or a licensed electrical heating technique, the effluent is passed through a heating zone, held at temperature and then passed back through a regeneration and a cooling zone. Only on confirmation of success of all of these steps is the effluent cleared for release.

Using robust, proven design principles, the systems deal with **CL Biowaste Level 1 to 3** and take into account two main areas of concern. Firstly, the systems effectively sterilise or inactivate any harmful pathogens in the waste stream and secondly total containment must be assured at all times.

#### Applications

Biologics Production  
Vaccine production  
Animal Production  
Laboratories  
Research Institutions  
Animal Laboratories  
Research Laboratories  
Pharmaceutical facilities



# EDSflow™ DATASHEET

VERSION 2.4



## Welcome

Since our foundation in 1961, Suncombe has pioneered the development of innovative solutions for cleaning in place, bio-waste decontamination, GMP Washers, sanitary skids and vessel skids. The business continues to be privately owned and managed day to day by Dave Adams and Steve Overton.

Supporting Dave and Steve is a close-knit, dedicated, highly motivated and long-standing team encompassing a wealth of technical experience and knowledge in all relevant disciplines, including design, manufacture, testing, installation, validation, documentation and after-sales support. All of our work is carried out across our own facilities, just off the M25 in north London.


The team employ the very latest techniques, standards and best in class solutions. Having such a strong team allows us to offer the ability to carry out all of our work in-house, under our direct control

and quality management system. It also ensures that we own and preserve all the knowledge and experience gained with every project and allows us to offer continued support for all our installed systems throughout their lifetime.

## Our Clientele



## Our Equipment



### Suncombe Ltd

Jade House, Lockfield Avenue, Brimsdown, Enfield, Middlesex, EN37JY, United Kingdom

T +44(0)20-8443-3454 F +44(0)20-8443-3969

E info@suncombe.com W www.suncombe.com

# EDSflow™ DATASHEET

VERSION 2.4

## Key Features

- Robust and reliable
- Pre-tested at Suncombe
- Plug n Play after delivery
- Steam utility required, integral steam generator or electrical heating
- Adjustable kill temperature and time
- Pumped discharge to sewers

Containment Level	BSL1, BSL2, BSL3
Treatment Vessels	None – continuous
Collection Vessels	Any Number
System Capacity (litres)	Up to 180,000L per day
Treatment Parameters	Variable f0 values and temperatures up to 140°C
System Cooling	>42C
Material of Construction	316 Stainless Steel
Operating Voltage	400 VAC
Automation	BioSuite Level 2000 standalone control system with touch panel display.



EDSflow™ Image

## Range

Part #	Capacity (litres per hour/Litres per day)
EDSflow300	300/5400
EDSflow500	500/9000
EDSflow1000	1000/18000
EDSflow2000	2000/36000
EDSflow3000	3000/54000
EDSflow5000	5000/90000
EDSflow10000	10000/180000

## Options

Material of Construction	SAF 2205 duplex or Hastelloy for chlorine resistance
Operating Voltage	Various
Automation	BioSuite Level 3000 with 21CFR 11 reports, records and networking
Safety	Failsafe versions available for SIL 2 and SIL 3
Covers	Thermoplastic or stainless steel covers
CIP	Chemical CIP for system cleaning
Pump Feed	Discharge pump to transfer waste
Transfer Method	Pumped, vacuum, eductor, differential pressure
Remote Control	Remote start and remote HMIs

Key Features	Benefits
Fully automated continuous processing	Provides ultra-efficient biowaste effluent Decontamination, enabling fully automated treatment of waste and safe release to drain.
Sustainable Design	The EDS incorporates facilities to recover up to 75% of the input energy.
Utilities	Steam utility required or integral steam generator or electrical heating + electrical supply, water and compressed air
Sanitary 316L stainless steel construction and components	Systems are constructed to the highest sanitary standards with 3.1/2.2 material traceability and welding dossier. This ensures a fully validatable and cleanable treatment environment.
Treatment Parameters	Variable configuration to provide Temperature/Time, f0 or log kill requirements. Temperatures up to 140°C.
Collection and Discharge Buffer Options	For large inlet flows, Integral Collection Tanks can be provided and for outlet flow limitations, integral discharge tanks can also be provided.
Vent Filter	Sterile, HEPA Vent Filter included, single, dual and duplex parallel arrangements available
Continuous monitoring of key parameters	Decontamination process is highly repeatable and validatable.
Containment Level	Systems suitable for the treatment of BSL 1, 2 and 3 Waste.
Safety	Alarms, interlocks and fail-safe design prevent waste discharge in the event of an unsuccessful treatment. This encompasses scenarios such as power loss and under-temperature events.
Siemens PLC and 12" colour HMI with options for larger HMIs	Control hardware is industry standard and supported worldwide by Siemens. Ethernet interface included for transfer of critical operating variables to other systems. Designed to enable integration to third party equipment or higher level control system.
Suncombe BioSuite software	Control software has been developed and proven over many years for EDS applications and includes a wide range of user or administrator configurable parameters to enable customised decontamination profiles. User passwords, Active Directory, Audit Trails, Electronic batch reports for local or network storage are possible. User interface screens and process visualisation is simple, intuitive, clear and comprehensive. Remote access options are possible if required. Software complies with FDA 21CFR and EU GMP regulations.
Fully automated report	Electronic pdf reporting included – printed report optional
Automatic Operation	Automatic Waste Inlet, treatment and Discharge Automatic Alarms and Warnings
Coolant Utility	Water or glycol Coolant.
Drainage Temperature	<42°C
Configurable	Based on standard modules, we can supply individual units custom designed for your specific requirement.
Cleaning In Place	Prepared for CIP with optional Automatic or Manual Cleaning In Place System
Waste Distribution	Optional sump pumps and discharge Stations
Redundancy	System configurations available for N+1 Dual redundancy

# EDSflow™ DATASHEET

VERSION 2.4



## EDSflow™ Configuration

The process flow diagram below depicts the typical EDSflow configuration and include; a) Single or dual start-up/CIP tanks b) One, two, three or four chemicals, c) Single or duty/standby transfer pumps and strainers d) Sump pump transfer from local sump. Using robust, proven design principles, the system uses a combination of thermal, elevated pressure processes designed to deactivate and sterilise waste, while ensuring absolute containment and protection of the external environment, ensuring that the target microbiological agents in the effluent are inactive before release.

## EDSflow™ Process Flow Diagram

